

FreeRider®

Luggie Elite Plus

Start your journey today with a limited edition Luggie Elite Plus. The deluxe model is a one of kind scooter with all the best features of the Luggie Standard and Elite. It comes with a wider seat, larger 10.5 lithium ion battery (for longer range), and comfortable armrests. It is lightweight, portable, and the best traveling scooter out on the market. Say "I can" and gain your independence back today!

Smart fold, Easy go

Foldable in 3 easy steps

Solid PU tires

Padded seating

Adjustable tiller

SPECIFICATION

Length : 39"/26"(folded)

Width : 18"

Weight(Scooter) : 51.7lbs

Weight(Scooter+Batteries) : 56.9lbs

Weight capacity : 320lbs

Turning radius : 37.5"

Ground clearance : 2.5"

Maximum speed : 4mph

Battery : 10.5AH Lithium Battery

Maximum incline : 6degrees

Wheels/Tires : Front Wheels:2"x6" PU solid

Rear Wheels:2"x8" PU solid

Drive system : 24-volt DC motor, Sealed

mini-transaxle, rear-wheel drive



COLOURS

Red

Blue

Black

